



DIAMOND KNIVES FOR MICRO-INCISION CATARACT SURGERY

SPECIFICALLY DESIGNED FOR MICS, DIAMATRIX'S GEM-QUALITY DIAMOND KNIVES CAN DELIVER A PRECISE INCISION TO ACHIEVE THE ARCHITECTURE YOU REQUIRE.

REALIZE MINIMAL INDUCED ASTIGMATISM ◆

INFLUENCE A DECREASE IN WOUND HEALING TIME ◆

ACHIEVE TIGHT SEALING INCISIONS FOR CHAMBER BALANCE ◆

PLUS THESE POPULAR STYLES:

1.0/ 1.3MM TAPERED B.U.C.K

2.0MM TRIAMOND

1.2/ 1.5MM MICROTRAP

COMPLIMENTARY EVALUATIONS ◆

CUSTOM SIZES AVAILABLE ◆

1544 Sawdust Rd, Suite 502
The Woodlands, TX USA 77380
www.diamatrix.com info@diamatrix.com
800-867-8081 (US) 281-367-8081 (T) 281-292-5481 (F)



Diamatrix Diamond Knives for MICS

Blade Configuration	Blade Width	Handle Type / Option			
		Premium Freehand		Stepped	
		Straight	Angled	Single	Multi
MicroTrapezoid	1.2/1.5mm	N 03880	N 04880		N 05880
MicroTrapezoid	1.5/1.8mm	N 03875	N 04875		N 05875
UltraTrapezoid	1.9/2.2mm	N 03860	N 04860		N 05860
Williamson Trapezoid	1.9/2.5mm	N 03870	N 04870	N 05870	
Clear Corneal	2.2mm	N 03600	N 04600	N 05600	
Clear Corneal	2.4mm	N 03605	N 04605	N 05605	
Shepherd MAC	2.2mm	N 03500	N 04500	N 05500	
Shepherd MAC	2.4mm	N 03505	N 04505	N 05505	
Tapered BUCK	1.0/1.3mm		N 04485		
Triamond	2.0mm		N 04450		